

Neutralon®

Neutralon® latex powdered surgical gloves are brown in color to reduce glare from surgical lighting and their micro-thin, polymer inner membrane allows for easy dry- and damp-hand donnability.

Safety & Protection

Infection Control

Viral Penetration. Although not an FDA requirement, Ansell's premium line of surgical gloves is tested to ensure protection against viral permeation. These gloves meet or exceed ASTM surgical glove standards. Manufactured within the quality guidelines of ISO 13485, ISO 9001, and FDA-QSR.

Cuff Roll-Down. This glove features an inverted, beaded cuff that prevents cuff roll-down.

Bone Cement. Always follow manufacturers' instructions when mixing and handling bone cement. This glove provides a protective barrier for handling and mixing bone-cement product. Double gloving is the usual practice.



Skin Protection

Low Protein Content

Lowry Protein Assay	Typical Results (µg/g)
Avg.	61
Min.	46
Max.	80

*This glove contains 60 µg/g or less of water/extractable protein per gram.
Caution: Safe use of these gloves by latex-sensitized individuals has not been established.*

Comfort & Fit

Donning

Alcohol Resistant. Ansell latex gloves provide a barrier for splash protection and/or short-duration of alcohol contact.

Damp-Hand Donnability. This glove contains a USP absorbable dusting powder that allows for smooth donning during dry or damp applications.

Double Donnability. This glove can be used as an under glove or an over glove depending on the user's preference.

Hand Fatigue

Shape of Glove/Former. Designed with Anatomic Curved Former for a more natural fit and less fatigue. This exclusive Ansell former is the result of a staff dedicated to the development and continuous improvement of formers.

Stretchability (Modulus). Because this glove is latex, wearers experience a more secure fit that provides greater comfort without being baggy.

Tactile Sensitivity

Dependability. This glove is recommended for times when you need more sensitivity to very low levels of tactile feedback.

Grip

Application Driven. Anti-slip, smooth finish for a comfortable, secure grip.

Caution: This product contains natural rubber latex, which may cause allergic reactions.

NEUTRALON®



Ansell

Everything you touch...we touch.