

Material Natural rubber from latex in compliance with ASTM D1076

Packaging 1 Sterile Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

Product Numbers

| SIZE (identified inside glove) | ASPN (Ansell Standard Product Number) |
|-----------------------------------|--|
| 5 1/2 | 5355 |
| 6 | 5360 |
| 6 1/2 | 5365 |
| 7 | 5370 |
| 7 1/2 | 5375 |
| 8 | 5380 |
| 8 1/2 | 5385 |
| 9 | 5390 |

Dimensional Specifications (minimum length 280mm)

| SIZE | Average Length mm (inches) | Average Palm Width mm (inches) |
|-------|-------------------------------|-----------------------------------|
| 5 1/2 | 287 (11.3) | 72 (2.8) |
| 6 | 287 (11.3) | 77 (3.0) |
| 6 1/2 | 287 (11.3) | 83 (3.3) |
| 7 | 287 (11.3) | 89 (3.5) |
| 7 1/2 | 287 (11.3) | 95 (3.7) |
| 8 | 287 (11.3) | 102 (4.0) |
| 8 1/2 | 287 (11.3) | 108 (4.3) |
| 9 | 287 (11.3) | 114 (4.5) |

Thickness (single)

| Data from Sampling | Average Thickness |
|----------------------|-------------------|
| Fingertip, mm (mils) | 0.26 (10.2) |
| Palm, mm (mils) | 0.23 (8.9) |
| Cuff, mm (mils) | 0.17 (6.5) |

Physical Properties (ASTM, EN, and ISO Requirements)

| | Before Aging | After Aging |
|---|--------------|-------------|
| Freedom from Holes (inspection level I) | 1.5 AQL* | N/A |
| Tensile, MPa (p.s.i.) minimum | 24 (3,480) | 18 (2,610) |
| Elongation (%) minimum | 750 | 560 |
| Modulus, MPa (p.s.i.) maximum | 5.5 (800) | N/A |

Ansell manufacturing specifications meet or exceed all the above.

| | |
|--------------------------------|---|
| Color | Brown |
| Accelerator | Activated aqueous dithiocarbamate |
| Donning Agent | Absorbable dusting powder, USP |
| Residual Powder Content | <15 mg/dm ² |
| Certification | Manufacturing facility registered to ISO 9001, ISO 13485, FDA-QSR, and CE certification |
| Quality Control | In compliance with ASTM D 3577, Type 1 |
| Sterilization | Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000 |
| References | ANSI/ASQC Z1.4 Sampling Procedures and Tables for Inspection by Attributes |
| Country of Origin | Malaysia |
| Storage | Cool environment, shielded from ozone and UV light |

* Acceptable Quality Level – maximum % of gloves with defects

Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.



For additional information and samples, call Ansell Healthcare at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
 200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
 Ansell and Neutralon are trademarks owned by Ansell Limited or one of its affiliates. ©2008 All Rights Reserved.