

Material Natural rubber from latex in compliance with ASTM D1076-06 containing dispersed bismuth oxide.

Packaging 1 Sterile Pair/Package – 5 Pairs/Case

Product Numbers

| SIZE (identified inside glove) | ASPN (Ansell Standard Product Number) |
|-----------------------------------|--|
| 6 1/2 | 20873065 |
| 7 | 20873070 |
| 7 1/2 | 20873075 |
| 8 | 20873080 |
| 8 1/2 | 20873085 |
| 9 | 20873090 |

Dimensional Specifications (minimum length 280mm)

| SIZE | Average Length mm (inches) | Average Palm Width mm (inches) |
|-------|-------------------------------|-----------------------------------|
| 6 1/2 | 280 (11.0) | 83 (3.3) |
| 7 | 280 (11.0) | 89 (3.5) |
| 7 1/2 | 285 (11.2) | 95 (3.7) |
| 8 | 285 (11.2) | 102 (4.0) |
| 8 1/2 | 290 (11.4) | 108 (4.3) |
| 9 | 290 (11.4) | 114 (4.5) |

Thickness (single)

| Data from Sampling | Average Thickness |
|----------------------|-------------------|
| Fingertip, mm (mils) | 0.35 (14) |
| Cuff, mm (mils) | 0.30 (12) |

Physical Properties (ASTM, EN, and ISO Requirements)

| | Before Aging | After Aging |
|---|--------------|-------------|
| Freedom from Holes (inspection level I) | 1.5 AQL* | N/A |
| Tensile, MPa (p.s.i.) minimum | 19 (2,755) | 17 (2,465) |
| Elongation (%) minimum | 730 | 560 |
| Force at Break (N) | 21 | 18 |

Ansell manufacturing specifications meet or exceed all the above.

| | |
|--------------------------------|---|
| Color | Khaki |
| Accelerators | Zinc diethyldithiocarbamate (ZDBC) Zinc dibutyldithiocarbamate (ZDEC) |
| Donning Agent | None (powder-free) |
| Residual Powder Content | < 2 mg/glove |
| Certification | Manufacturing facility registered to ISO 9001 and ISO 13485 Contains 50 micrograms or less of total water-extractable protein per gram per glove |
| Quality Control | In compliance with ASTM D3578-00 |
| Sterilization | Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000 |
| References | ISO 2859-1 Sampling Procedures and Tables for Inspection by Attributes |
| Country of Origin | Malaysia |
| Storage | Cool environment, shielded from ozone and UV light |

*Acceptable Quality Level – maximum % of gloves with defects

Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.



For additional information and samples, call Ansell Healthcare at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell is a registered trademark owned by Ansell Limited or one of its affiliates. ©2007 All Rights Reserved.

SS390 11/2007

Ansell
Everything you touch...we touch.