



Give your team
more control.



INTRODUCING

Biogel® PI UltraTouch™ G
Biogel® PI UltraTouch™ M with textured surface

NEW GLOVES. INCREASED GRIP.

Two high-quality additions to our complete line of NON-LATEX gloves.



SETTING THE STANDARD IN PROTECTION

Get MORE GRIP CHOICES as two new NON-LATEX surgical gloves join the Biogel® polyisoprene product line.

INTRODUCING Biogel® PI UltraTouch™ G and Biogel® PI UltraTouch™ M.

Both gloves offer the fit, feel and comfort of latex while providing optimal grip for a range of surgical procedures.

Biogel PI UltraTouch G

- Polyisoprene
- Micro-roughened surface
- Curved-finger design



Biogel PI UltraTouch M

- Polyisoprene
- Textured surface
- Straight-finger design



ACTUAL COLOR

RE-ORDER # 421

#	SIZE	PAIRS
42155	5 ½	50/Box
42160	6	50/Box
42165	6 ½	50/Box
42170	7	50/Box

4 BOXES PER CASE



ACTUAL COLOR

RE-ORDER # 426

#	SIZE	PAIRS
42655	5 ½	50/Box
42660	6	50/Box
42665	6 ½	50/Box
42670	7	50/Box

4 BOXES PER CASE

#	SIZE	PAIRS
42675	7 ½	50/Box
42680	8	50/Box
42685	8 ½	50/Box
42690	9	40/Box

Both gloves are recommended for use in microvascular, cardiothoracic, ENT, plastic, neurosurgical, gynecological and general surgical procedures or surgeries where latex allergies are a concern for patients or clinicians.

WHY CHOOSE BIOGEL?

- Every glove is air inflation-tested and visually inspected for quality and user confidence
- Low dermatitis potential
- Industry-leading AQL freedom from holes of 0.65¹
- Competitors' gloves are typically twice as likely to have a hole²
- Biogel is the only major surgical glove brand with a non-pyrogenic range³
- Biogel gloves exceed industry standards for⁴
 - Force at break
 - Tensile strength
 - Elongation

Biogel PI UltraTouch M and Biogel PI UltraTouch G can be used alone or with the Biogel® PI Indicator™ Underglove to create the Patented Puncture Indication System™.



SETTING THE STANDARD IN PROTECTION



The Mölnlycke Health Care and Biogel® names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Distributed by Mölnlycke Health Care US, LLC, Norcross, Georgia 30092. ©2010 Mölnlycke Health Care AB. All rights reserved.

1 0.65 AQL for Freedom from Holes ASTM D3577, Test Method D5151

2 "Implementation recommendations for making health care facilities latex safe." Kim K.T., Graves P.B., Safadi G.S., Alhadeff G., Metcalfe J., AORN Journal, March 1998

3 Dorland W.A. Dorland's Illustrated Medical Dictionary, 30th Edition. Philadelphia: Saunders, 2003.

4 MKT004 Why Choose Biogel