



Biogel® Eclipse®


A softer more supple glove than standard latex gloves



- Smooth surface
- Outstanding fit, feel and comfort, with a softer feel than standard latex gloves.
- Biogel coating for ease of donning with damp or dry hands
- 0.65 AQL for holes¹
- 100% air inflation and visually inspected²
- Meet or exceed ASTM surgical glove standards
- Beaded cuff
- Sure-Flap packaging

Product Properties

- Material: Natural rubber latex with Biogel polymer inner coating
- Biogel coating soothes the skin to help prevent moisture loss
- Nominal thickness: fingertip 0.22mm (9mils)
- Curved finger design
- Non-pyrogenic
- Low dermatitis potential³
- Contains 50 µg or less of water extractable protein per gram⁴
- Color: Straw



| Order # | Size | Pairs/Box |
|---------|------|-----------|
| 75255 | 5½ | 50 |
| 75260 | 6 | 50 |
| 75265 | 6½ | 50 |
| 75270 | 7 | 50 |
| 75275 | 7½ | 50 |
| 75280 | 8 | 50 |
| 75285 | 8½ | 50 |
| 75290 | 9 | 40 |

4 Boxes per Case

¹0.65 AQL for freedom from holes. ASTM D3577, Test Method D5151.

² Each Biogel glove undergoes multiple checks for quality and user confidence.

³ Low dermatitis – gloves demonstrated reduced potential for sensitizing user to chemical additives.

⁴ This product contains 50 µg or less of water extractable protein per gram. **CAUTION: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established. WARNING: Do not use this product if you have a known allergy to natural rubber protein or chemical additives.**