

...“Double gloving with an **indicator system provides the best protection** and allows the timeliest identification of perforations!”

Only Biogel's double gloving **Patented* Puncture Indication System** clearly alerts a practitioner if and where glove perforations occur in the presence of fluid, so you can just focus on the surgery at hand.

The Biogel double gloving **Patented* Puncture Indication System** has a 97% accuracy rate in identifying punctures.²

More than 90% of the participants in a recent study rated the **Patented* Puncture Indication System** good or better for fit, comfort, sensitivity, and speed of indication than other single and double gloves tested.³

Works with a variety of Biogel translucent surgical gloves, which allows you to create your own customized double gloving **Patented* Puncture Indication System**.

“Reported glove perforation rates are as high as 78% in high risk procedures, and even laparoscopic procedures report around 20%. Double gloving significantly decreases exposure to blood and provides better protection.”¹

The American College of Surgeons Committee on Operating Room Environment claims, “Double gloving does help to cut down by a factor of ten (10) the number of potential exposures...”⁵

Spot the risk. Call 1.800.516.9236 for a FREE sample.

*U.S. Patent Nos. 5,224,221 and 5,524,294



Biogel® Patented Puncture Indication System

spot the risk.

NON-LATEX
Now Available



MÖLNLYCKE®
HEALTH CARE

¹Florman S et. al. Efficacy of Double Gloving with an Intrinsic Indicator System. Surgical Infections 2005; 6(4):385-395. ²Wigmore SJ and Rainey JB. Use of Coloured Undergloves to Detect Glove Puncture. BJS 1994; 81:1480
³Florman S et. al. Efficacy of Double Gloving with an Intrinsic Indicator System. Surgical Infections 2005; 6(4):385-395. – Biogel Indicator Underglove with Biogel Super-Sensitive as over glove ⁴ACS. Preventing Disease in the Operating Room. Panel discussion, American College of Surgeons Spring Meeting, April 29, 1998. Accessed at http://www.facs.org:80/about_college/acscomm/core/coreprev.html on September 16, 1999.

Biogel, the Biogel logo and Mölnlycke are registered trademarks of Mölnlycke Group of companies. Distributed by Mölnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2006 Regent Medical Limited. All rights reserved. LIT# 1-0141

ONLY Biogel gloves paired together can offer the exclusive double gloving **Patented Puncture Indication System** – available in a variety of latex and non-latex, Biogel has the customized fit and comfort that is just right for you.

Undergloves

Biogel® Indicator™ Underglove

the patented puncture indication system, designed to alert wearer if a perforation occurs in the presence of fluid



Order #	Size	Pairs/Box*
31255	5½	50
31260	6	50
31265	6½	50
31270	7	50
31275	7½	50
31280	8	50
31285	8½	50
31290	9	40

*4 boxes per case.

Outergloves

Biogel® Eclipse™

made with deproteinised latex, a softer more supple gloving category than standard latex gloves



Order #	Size	Pairs/Box*
40655	5½	50
40660	6	50
40665	6½	50
40670	7	50
40675	7½	50
40680	8	50
40685	8½	50
40690	9	40

*4 boxes per case.

Biogel®

the world's finest general use surgical glove



Order #	Size	Pairs/Box*
30455	5½	50
30460	6	50
30465	6½	50
30470	7	50
30475	7½	50
30480	8	50
30485	8½	50
30490	9	40

*4 boxes per case.

Biogel® Optifit™ Orthopaedic

approximately 20% thicker-than Biogel with ultra-precise anatomically engineered fit gloves



Order #	Size	Pairs/Box*
31060	6	40
31065	6½	40
31070	7	40
31075	7½	40
31080	8	40
31085	8½	40
31090	9	40

*4 boxes per case.

Biogel® M

same as Biogel surgical glove but with micro-textured surface for enhanced grip and wider palm



Order #	Size	Pairs/Box*
30555	5½	50
30560	6	50
30565	6½	50
30570	7	50
30575	7½	50
30580	8	50
30585	8½	50
30590	9	40

*4 boxes per case.

Other latex Biogel outer glove options include: Biogel Sensor®, Biogel® Super-Sensitive™, and Biogel® Optifit™ Super-Sensitive™.

NON-LATEX

Biogel® Skinsense™ Underglove

can be used alone or as a non-latex underglove in patented puncture indication system



Order #	Size	Pairs/Box*
40655	5½	50
40660	6	50
40665	6½	50
40670	7	50
40675	7½	50
40680	8	50
40685	8½	50
40690	9	40

*4 boxes per case.

Biogel® Skinsense™

Synthetic surgical glove that performs like latex.



Order #	Size	Pairs/Box*
31455	5½	50
31460	6	50
31465	6½	50
31470	7	50
31475	7½	50
31480	8	50
31485	8½	50
31490	9	40

*4 boxes per case.

Double gloving

with the Patented* Puncture Indication System



1. Choose from our latex or non-latex underglove.
2. Then don your preferred Biogel outerglove.

3. A dark spot enables rapid detection of glove perforation in the presence of fluid.



Biogel® Patented Puncture Indication System

spot the risk.



*U.S. Patent Nos. 5,224,221 and 5,524,294

Biogel, the Biogel logo and Mölnlycke are registered trademarks of Mölnlycke Group of companies. Distributed by Mölnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2006 Regent Medical Limited. All rights reserved. LIT# 1-0141