



Biogel® PI



The synthetic polyisoprene alternative that's hard to believe it's not natural latex



- Synthetic polyisoprene material shares similar molecular structure with natural rubber latex
- Can be used alone or as the outer glove, in combination with Biogel® SkinSense® Underglove, for our latex-free **Patented Puncture Indication System**
- Powder-free to eliminate starch powder-related complications
- Biogel coating for ease of donning with damp or dry hands
- 0.65 AQL for holes¹
- 100% air inflation and visually inspected²
- Meet or exceed ASTM surgical glove standards
- Beaded cuff
- Sure-Flap packaging
- Non-pyrogenic

Product Properties

- Material: Synthetic polyisoprene with Biogel polymer inner coating
- Biogel coating soothes the skin to help prevent moisture loss
- Nominal thickness: fingertip 0.30mm (12 mils)
- Non-pyrogenic
- Low dermatitis potential³
- Curved finger design
- Color: Cream



Order #	Size	Pairs/Box
40855	5½	50
40860	6	50
40865	6½	50
40870	7	50
40875	7½	50
40880	8	50
40885	8½	50
40890	9	40

4 Boxes per Case

¹0.65 AQL for freedom from holes ASTM D3577, Test Method D5151.

²Each Biogel glove undergoes multiple checks for quality and user confidence.

³Low dermatitis potential – gloves demonstrated reduced potential for sensitizing user to chemical additives.