



Biogel Sensor[®]


Thinner for improved sensitivity



- Approximately 20% thinner than Biogel M
- Recommended for procedures where sensitivity and fingertip control are crucial
- Specially formulated micro-textured surface for positive grip
- Powder-free to eliminate starch powder-related complications
- Biogel coating for ease of donning with damp or dry hands
- 0.65 AQL for holes¹
- 100% air inflation and visually inspected²
- Meet or exceed ASTM surgical glove standards
- Beaded cuff
- Sure-Flap packaging

Product Properties

- Material: Natural Rubber Latex with Biogel polymer inner coating
- Biogel coating soothes the skin to help prevent moisture loss
- Nominal thickness: fingertip 0.21mm (8mils)
- Straight finger design
- Non-pyrogenic
- Low dermatitis potential³
- Contains 50 µg or less of water extractable protein per gram⁴
- Color: Straw



Order #	Size	Pairs/Box
30655	5½	50
30660	6	50
30665	6½	50
30670	7	50
30675	7½	50
30680	8	50
30685	8½	50
30690	9	40

4 Boxes per Case

¹ Known allergy to natural rubber protein or chemical additives.
² No reactions. Safe use of this product by or on latex sensitized individuals has not been established.