

# Biogel® Skinsense® Underglove

Can be used alone or as a non-latex underglove




- Smooth surface for use alone or as an underglove with any Biogel non-latex outer glove combination
- Chloroprene composition eliminates possibility of glove-related latex protein sensitization
- Powder-free to eliminate starch powder-related complications
- Biogel coating for ease of donning with damp or dry hands
- 0.65 AQL for holes<sup>1</sup>
- 100% air inflation and visually inspected<sup>2</sup>
- Meet or exceed ASTM surgical glove standards
- Beaded cuff
- Sure-Flap packaging

## Product Properties

- Material: Chloroprene with Biogel polymer inner coating
- Biogel coating soothes the skin to help prevent moisture loss
- Nominal thickness: fingertip 0.20mm (8mils)
- Curved finger design
- Low dermatitis potential<sup>3</sup>
- Color: Dark Blue tone

Order #	Size	Pairs/Box
40655	5½	50
40660	6	50
40665	6½	50
40670	7	50
40675	7½	50
40680	8	50
40685	8½	50
40690	9	40

4 Boxes per Case



<sup>1</sup>0.65 AQL for freedom from holes ASTM D3577, Test Method D5151.

<sup>2</sup>Each Biogel glove undergoes multiple checks for quality and user confidence.

<sup>3</sup>Low dermatitis potential – gloves demonstrated reduced potential for sensitizing user to chemical additives.