



DURAPRENE

Synthetic
freedom

Latex-like
performance

Outstanding
durability

Comfortable fit

Powdered or
powder-free

Protection, Performance and Peace of Mind.

 **Allegiance**

The Synthetic Glove for Authentic Peace of Mind.



If you prefer a latex-free glove, there's no reason to compromise protection and performance. Our Duraprene™ Synthetic Surgical Glove is produced with an innovative synthetic formulation that fits, feels and performs like natural rubber latex. This advanced glove ensures strong confident barrier protection while offering outstanding comfort and tactile sensitivity.

Duraprene is available with or without powder. And with our patented natural anatomical design, this glove provides comfort and helps reduce hand fatigue in even the longest procedures.

Duraprene – get protection, performance and peace of mind from one synthetic surgical glove.

Ordering Information

Powder-Free Duraprene

Catalog Number	Size
2D72B01	5½
2D72B11	6
2D72B21	6½
2D72B31	7
2D72B41	7½
2D72B51	8
2D72B61	8½
2D72B71	9

Duraprene

Catalog Number	Size
2D72A01	5½
2D72A11	6
2D72A21	6½
2D72A31	7
2D72A41	7½
2D72A51	8
2D72A61	8½
2D72A71	9

All products packaged
40 per box/5 boxes per case

- Available powdered and powder-free.
- Patented synthetic formulation emulates natural rubber latex.
- Provides superior durability and barrier protection.
- Contains no natural rubber latex.
- Fit and feel of natural rubber latex in a latex-free glove.
- Firm, confident grip with an anti-slip finish and excellent tactile sensitivity.
- Meets or exceeds ASTM D3577 strength requirements.
- Beaded 12" cuff prevents rolldown and provides extended hand and wrist coverage.
- Manufactured in a facility certified to ISO 9001 standards.

Professional Services at Allegiance

Our Professional Services Department offers ongoing service and support to ensure customer satisfaction. Our toll-free number provides a direct source for inquiries regarding clinical application and product usage. The Professional Services Group is staffed by registered nurses who offer a full range of clinical support services and product information.

To contact the Professional Services Department, call 1-800-327-7503.

To order, call Allegiance Customer
Service at 1-800-964-5227.



Allegiance Healthcare Corporation
Gloves
1500 Waukegan Road
McGaw Park, Illinois 60085
www.allegiance.net