

## G-VIR® > Description

<b>GENERAL INFORMATIONS</b>	
<b>Product</b>	Synthetic, powder-free, sterile surgical glove incorporating a disinfecting liquid, reducing the transmitted viral load in the event of a Blood Exposure Accident
<b>Packaging</b>	20 pairs per box
	3 boxes per carton
<b>Size</b>	6, 6 ½, 7, 7 ½, 8, 8 ½, 9
<b>SPÉCIFICATIONS</b>	
<b>Glove length</b>	260 mm minimum
<b>Glove thickness</b>	500±50µm
<b>Basic material</b>	Thermoplastic elastomer - <b>Natural latex-free</b>
<b>Biocompatibility</b>	Evaluation according to the ISO 10993
<b>Donning Layer</b>	<b>Powder-free</b> , Conforms to NF EN 455-3
<b>Sterilization</b>	Gamma radiation
<b>Use-by date</b>	4 years from the manufacturing date
<b>CERTIFICATIONS</b>	
<b>CE Mark</b>	Certified since June 2003
<b>ISO 13485:2003</b>	Certified since June 2003
<b>ISO 9001:2000</b>	Certified since November 2005
<b>Notified Body</b>	SGS United Kingdom Ltd Systems & Services Certifications (0120)
<b>MANUFACTURER'S DETAILS</b>	
<b>Address</b>	HUTCHINSON SANTE
	<b>Consumer Service</b> : Rue Marret et Paturel - 60140 Liancourt - France
	<b>Headquarters</b> : 2 rue Balzac - 75008 Paris - France
<b>Tel</b>	From France : 0820 167 560
	From abroad : + 33 1 72 87 48 08
<b>E-Mail</b>	<a href="mailto:info.gvir@gvir.com">info.gvir@gvir.com</a>